

**IMPRESSION
MATERIAL**



**CROWN &
BRIDGE
MATERIAL**

**RESTORATIVE
AND
FILLING**



**PRODUCT
CATALOGUE
2025-26**



**BONDING
ADHESIVE
& DESENSITIZER**

**CEMENTS
AND
LINERS**



The *ALTERNATIVE*
Premium Brand
From the **U.K.**

innovation with passion

'complete commitment to quality'



Since 2004, **Medicept UK Ltd** has been delivering products that ensure excellence in dental care, and become the trusted choice of satisfied dentists worldwide:

Albania • Algeria • Azerbaijan • Bahrain • Bangladesh • Bulgaria
Burkina Faso • Colombia • Costa Rica • Cyprus • Egypt • Ethiopia
Greece • Honduras • India • Indonesia • Ireland • Israel • Italy
Kenya • Kosovo • Kuwait • Libya • Lithuania • Malaysia • Malta
Moldova • Macedonia • Nepal • Pakistan • Philippines • Poland
Portugal • Qatar • Romania • Saudi Arabia • Sri Lanka • Spain
Tanzania • Turkey • United Kingdom • Uruguay • Uzbekistan • Yemen

XtraLute

GLASS IONOMER LUTING CEMENT



- Exceptional strength - Long lasting cementation
- Good radiopacity - Shows clearly under x-ray
- Easy mixing - To reduce surgery time
- Excellent biocompatibility - Minimal pupal reaction
- Compressive Strength 122MPa

Ref No : 2003 - 15g Powder / 10ml Liquid
Ref No : 2004 - 35g Powder / 20ml Liquid

XtraCem

GLASS IONOMER RESTORATIVE CEMENT



- Excellent early translucency - Superior aesthetics
- Exceptional strength - Long lasting restoration
- Packable consistency - For ease of handling
- Easy mixing - To reduce surgery time
- Compressive Strength 228 MPa

Ref No : 2001 - 15g Powder / 10ml Liquid

XtraCem-P

RADIOPAQUE GLASS IONOMER POSTERIOR RESTORATIVE CEMENT



- Superb translucency - Excellent appearance
- Aesthetic match to natural tooth - Visually superior
- High compressive strength - Durable restorations
- Adheres chemically to tooth structure - No need to etch or bond
- Compressive Strength 260 MPa.

Ref No : 2006 - 15g Powder / 10 ml Liquid

XtraLuteRMGI

RESIN MODIFIED GLASS IONOMER LUTING CEMENT



This resin modified glass ionomer luting cement is ideal for luting crowns, bridges, inlays, onlays and "all ceramic" restorations.

- High strength for lasting cementation
- Colour change makes trimming easy
- Aesthetic universal shade

Ref No : 2005 - Powder 15gm / Liquid 7ml, mixing pad & scoop

XtraCem-LC

RADIOPAQUE LIGHT-CURED GLASS IONOMER CEMENT



A light curing, resin reinforced glass ionomer restorative. The material features both chemical and light curing mechanisms to ensure the best mix of early strength and command setting. Offering both excellent compressive and flexural strengths to provide durable restorations.

- Compressive Strength 221MPa

Ref No : 2002 - 12g Powder / 5ml Liquid

ZinCem

ZINC PHOSPHATE CEMENT



- Adhesive to enamel, dentine and non precious metals
- Wide range of uses
- Proven formula - Reliability assured
- Radiopaque - Shows clearly under x-ray
- Easy to mix - To reduce surgery time
- Compressive Strength 148 Mpa

Ref No : 2008 - 90g Powder / 30 ml Liquid
Ref No : 2007 - 30g Powder / 15 ml Liquid

PolyCem

POLYCARBOXYLATE CEMENT



- Sustained fluoride release
- Prevention of secondary caries
- Easy mixing - To reduce surgery time
- Adhesive to enamel, dentine and non precious metals
- Wide range of uses
- Contains no phosphoric acid - Minimal pupal reaction
- Radiopaque - Shows clearly under x-ray
- Compressive Strength 96 Mpa

Ref No : 2010 - 90g Powder / 30 ml Liquid
Ref No : 2009 - 30g Powder / 15ml Liquid

HydroCal

CALCIUM HYDROXIDE PASTE



HYDROCAL is a Calcium Hydroxide based two paste system used as a liner in deep cavities and for direct and indirect pulp capping

- Convenient paste mixing product in tube/tube presentation
- Rigid, self-setting radiopaque composition
- Will not inhibit the polymerisation of acrylic and composite restoration
- Withstands condensation of amalgam

Ref No : 9005 - 13g Base Paste, 11g Catalyst Paste

HydroCal

LIGHT CURE CALCIUM HYDROXIDE CAVITY LINER



Hydrocal LC is a unique light-cured, radiopaque, dentinadhering liner and base material containing calcium hydroxide and calcium hydroxyapatite in a urethane dimethacrylate (UDMA) base. It is a highly filled resin (for a liner) with minimal shrinkage and water absorption. Hydrocal LC has a higher viscosity than conventional liners, eliminating the need for a ball-type applicator. Its superior strength and water insolubility make Hydrocal LC the material of choice under composite, amalgams & gold

Ref No : 5007 - 5g Syringe & 3tips

implalute IMPLANT CEMENT SEMI PERMANENT IMPLANT CEMENT



Implalute cement for superstructures on a resin base 4:1. Specially designed for semipermanent fixing on customised abutments, abutments with reduced adhesion surface or particularly small abutments. With reinforced adhesion and compressive strength. The restoration can still be removed without damage even after extended time in position. Compatible with all material combinations, dual curing system. Colour: natural opaque.

Ref No : 9002 - 5ml mini-mix Syringe, 5 Mixing Tips

selfcem DUAL-CURE SELF-ADHESIVE RESIN CEMENT



Exceptional bond strength - one of the most retentive self-adhesive auto-mix cements available • Dual-cure-Light-cure or self-cure for maximum versatility • Convenient auto-mix delivery system - direct and accurate one-step placement with minimal waste - no trituration required! • Contains fluoride - provides a sustained fluoride release to minimize the potential for caries • Multiple applications - Perfect for PFM, Zirconia and Alumina-based restorations, and posts • Versatile shade selection - Choose from A2 and Translucent

Ref No : 5009 - 5ml Automix Dual-Barrel Syringe, 5 Tips
Shade : A2 & Translucent

Glass Ionomer Universal Liquid



A Universal Glass Ionomer Liquid which can be used with glass ionomer restorative and luting.

Ref No : 2011 - 15ml Liquid



DUAL CURE COMPOSITE WITH ZIRCONIUM DIOXIDE FOR CORE BUILD-UPS AND POST CEMENTATION



- One material, two indications: core build-up and post-cementing
- Contains zirconium dioxide
- Very high flexural and compressive strength
- Excellent flow properties
- Elasticity behavior comparable to natural tooth substance
- AUTOMIX-Syringe

Ref No : 9013 - 2 x 5ml/9g Automix Syringes, 10 Mixing tips
Shade : Shade A3



NANO HYBRID COMPOSITE



Medicet's new restorative is an advanced composite for both anterior and posterior use.

- 2018** - Flash Starter Kit
3x4gm syg. 5ml Bond, 2g Etch, 50 microbrush
Shade: A1, A2, A3
- 3004** - Flash Syringe
Refill Pack 1 x 4g
Shade: A1, A2, A3, A3.5, B1, B2, OA2, OA3, B/W, Inc
- 3001** - Flash Kit
6 x 4 g Flash Syringes • 1 x 2 g SoloEtch Syringe
50 x Micro Applicators • 1 x 5 ml Bond Plus NT
1 x 2 g Flow Plus Syringe • Mixing Pad • Mixing Wells
Micro Tips for SoloEtch • Instruction Sheet
Shade: A1, A2, A3, A3.5, B1, OA2



FLOWABLE LIGHT CURE COMPOSITE



An advanced light curing flowable composite material. Offering controlled light diffusion for excellent aesthetics and visual appearance. Exceptional handling based on flow properties combined with superb strength and high elasticity for durable restorations. Radiopaque, unlike most other flowable composites.

Ref No : 3005 - 1 x 2g Syringe & 5 Tips
Shade: A1, A2, A3, A3.5, B2



TEMPORARY CAVITY FILLING MATERIAL



Ready-to-use filling material for temporary filling and sealing. Setting in the mouth under the effect of saliva, great adhesive strength in the cavity, excellent marginal sealing, easy to remove.
Colour: white opaque

Ref No : 5003 - 30g Jar



PIT & FISSURE SEALANT



- Thixotropic, Low-Viscosity - Able to flow into occlusal grooves
- 58% Filled, Low polymerization shrinkage, Superior strength
- Fluoride Releasing - Anti-caries benefit from sustained, fluoride release along margins of restoration
- Radiopaque - Visible on x-ray

Ref No : 5017 - 1.2ml Syringes, 3 Delivery Tips
Shade : Milk White Colour



PROPHYLAXIS PASTE



A blend of abrasives selected to ensure optimum removal of plaque and tartar without any appreciable abrasion to enamel. The thixotropic properties of the paste are carefully controlled to ensure easy application without undesirable spattering at the rotation speeds required for prophylaxis. Protec contains only natural flavourings, are fluoride free and are available in medium and course grades.

Ref No - 250g Jar - 6001 Orange, 6002 Mint, 6003 Strawberry, 6004 Bubblegum
50g Jar - 6005 Orange, 6006 Mint, 6007 Strawberry, 6008 Bubblegum

BOND *Plus* SE

SINGLE COMPONENT,
SELF-ETCH DENTAL ADHESIVE



A one step self etching bonding agent. No need to use etchants or primers prior to applying the adhesive. Just apply adhesive, light cure and you are ready to apply the composite restorative material. Bond strength and biocompatibility are both excellent and also the material is completely HEMA free. A true seventh generation bond.

Ref No : 4003 - 5ml Bottle

BOND *Plus* NT

TOTAL-ETCH SINGLE-COMPONENT
ADHESIVE SYSTEM



A single component prime and bond material, suitable for bonding composites and compomers, to enamel dentine and metal. The material is in the form of a viscous, non-thixotropic light curable liquid that has the ability to moisten and penetrate dentine.

Ref No : 4007 - 5ml Bottle

SoloEze

DESENSITIZER



SoloEze is a Hema Free Desensitizing solution. Prevents and rapidly alleviates enamel and dentinal hypersensitivity when filling and cementing, for temporary restorations, sensitive cervical areas.

Ref No : 5002 - 7ml bottle

SoloEtch

ETCHING GEL



Medicept Dental's Solo Etch 37% phosphoric is a medium viscosity etchant gel used for etching enamel and dentin prior to bonding procedures. Solo Etch is thixotropic, stays where placed and will not run. Does not dry out and has a long shelf life when properly stored. Washes off clean.

Ref No : 5005 - 50 ml Syringe & 50 Tips
Ref No : 5004 - 1 x 2g Syringe & 3 Tips

AccuCem^{Plus}

EUGENOL FREE TEMPORARY CEMENT



Totally Eugenol free cement based in zinc oxide and EBA. It gives a smooth and creamy consistency with a long working time.

- Radiopaque with low film thickness
- Does not inhibit resin polymerisation
- Ideal for temporary cementation of provisional crowns and bridges
- Suitable for use on patients with Eugenol sensitivity or allergy

Ref No : 9007P - 1 x 40g Base Paste & 1 x 20g Catalyst Paste

AccuTemp

TEMPORARY CROWN AND BRIDGE MATERIAL



Accutemp is a Bis-acrylic composite 10:1, self curing, for high-quality short- and long-term provisionals. Natural aesthetics and fluorescence, in all lighting conditions. Exceptional bending strength, optimal adaptation and fitting. Easy surface finishing (grinding, polishing). Tasteless and long-term colour stable during the whole use. Direct application with Automix2 cartridge system. Shade: A1, A2, A3.

Ref No : 9007 - 50ml Cartridge, 6 Mixing Tips
Shade : A1, A2, A3

Medi-Pak

KNITTED RETRACTION CORD



- Interlocking Knitted Loops - Unsurpassed retraction/displacement, Easy to control and place, Will not become entangled in burs or packers
- Non-impregnated, Holds more hemostatic solution than other cords
- Various Sizes - Useful for a variety of procedures and teeth.

Ref No : 5018 - (1) Bottle containing 8 feet (244 cm) of cord.
Available in the following sizes:
#000, #00, #0, #1, #2 & #3

HemoStat Gel

HEMOSTATIC GEL



- 20% ferric sulfate hemostatic gel
- Stops bleeding in seconds, saving chair time and decreasing the need for impression remakes.
- Hemostatic agent for crown and bridge procedures
- Viscous - Stays in place when delivered
- Unique Binding & Coating Agents
- Easy to Deliver with Infuser tip

Ref No : 5015 - 5ml Syringe, 5 Delivery Tips

HemoStat Clear

CLEAR HEMOSTATIC GEL



- 25% aluminum chloride formula rapidly eliminates minor bleeding.
- Recommended for anterior restorations,
- For use in the esthetic zone to prevent coagulum and hemostatic residue from adhering to restorations
- Colorless - Will not leave a colored residue after procedures in esthetic zone
- Viscous - Non-dripping

Ref No : 5016 - 5ml Syringes, 5 Delivery Tips

RC-Plus



ROOT CANAL PREPARATION CREAM FOR CHEMO-MECHANICAL PREPARATION OF ROOT CANALS.

RC-Plus unique formulation of glycol, urea peroxide and EDTA in a special water soluble base helps remove calcifications and lubricates the canal to permit more efficient instrumentation. Whether you prefer finger or engine driven instruments, RC-Plus allows reamers and files to move easily without binding. RC-Plus is also excellent for use with apex locators as it permits consistently reliable readings.

Ref No : 4002 - 9g Syringe

SureSeal



ZINC OXIDE EUGENOL-BASED ENDODONTIC CEMENT

- Barium Sulfate gives radiopacity for follow up.
- Sealing endodontic obturations.
- Cementation of Gutta Percha points or silver points, in the obturation of these root canals.
- Expands upon setting to create a positive seal.
- Powder/liquid presentation to adjust the cement viscosity to the technique required

Ref No : 3003 - 15ml bottle, 40gm powder
Ref No : 3021 - 10ml bottle, 15gm powder

SureSealXT



ROOT CANAL SEALER

Sure Seal XT is a paste/paste zinc oxide eugenol root canal sealer. It is light in color, non darkening, and radiopaque.

- Radiopaque • Barium sulfate • Low film thickness • Z O E - based formula • Low leakage and low dimensional change • Provides radiopacity while eliminating staining the tooth • Pastes mix easily to a smooth and creamy consistency for ease of use • For permanent obturation of the root canal space with the aid of obturation points

Ref No : 9004 - 8g Base Paste, 3.5g Accelerator Paste

selfpost
GLASS FIBER POST



Ref No :
8005 Pack of 10 Post - Size 0,1,2, 3 & 4
8006 Pack of 1 Drill - Size 0,1,2, 3 & 4
8007 Selfpost Starter Kit
3 Post each of size 0,1,2
1 Drill each of size 0,1,2

bio fill
Root Canal Filling Material



ROOT CANAL FILLING MATERIAL SILVER FREE

Biofill powder and resin are mixed to produce a root canal filling material of excellent sealing capacity. As a result of its close adaptation to the root canal walls and its very small contraction during setting, Biofill allows tightly sealing root canal obturation of excellent tissue compatibility.

Ref No : 4004 - 8g Powder/10g Resin

bio fill
Epoxy Resin Root Canal Sealer
silverfree



EPOXY RESIN ROOT CANAL SEALER PASTE

a base and resin mixed to produce a root canal filling material of excellent sealing capacity as a result of its close adaptation to the root canal walls and its small contraction during setting. Permanent obturation of root canals of teeth of the secondary dentition with the aid of obturation points.

- Excellent sealing capacity
- Easy to mix
- Paste/paste presentation

Ref No : 4004P - 20g Base Paste, 10g Catalyst Paste

Mixing Pad



Pack of 4

Ref No : 8003 - Pack of 4

Wave X



DENTAL FILM FOR INTRAORAL RADIOGRAPHY SPEED - E

Biofill powder and resin are mixed to produce a root canal filling material of excellent sealing capacity. As a result of its close adaptation to the root canal walls and its very small contraction during setting, Biofill allows tightly sealing root canal obturation of excellent tissue compatibility.

Ref No : 9014 - Pack of 150 Films

EXTREME

A-SILICONE IMPRESSION MATERIAL



- Two-Phase, one-step A-silicone impression system.
- Short intraoral setting time (2min)
- High consistency putty, regular setting.
- Low consistency wash, light bodied
- Super hydrophilic, thixotropic, crisp detail

Ref no: 5010 -Extreme Putty - 250ml base + 250ml catalyst
Ref no: 5011 - Extreme Lite - 2 x 50ml Cartridge
+ 8 Mixing tips

ORANGEBITE

BITE REGISTRATION MATERIAL



- A-silicone based, medium consistency - medium bodied, thixotropic.
- Extremely short time in the mouth (60 sec.),
- Maximum final hardness, approx. 45 Shore D, easy processing.
- Colour : Base orange, Catalyst white 1:1
- Flavour : Orange

Ref No : 5013 - 2 x 50ml Cartridge 8 Mixing Tips,
8 Bite Registration Tips.

EXTREME
MONO

MONO PHASE IMPRESSION SILICONE



Extreme Mono is a medium tray consistency for single step impression-taking or mixed technique, with distinct snap set.

Ref No : 9012 - 2 x 50ml Cartridge & 8 Mixing Tips

EXTREME
ADHESIVE

ADHESIVE FOR A SILICONE



For Adhesion of A-silicon impression material to plastic and custom trays, fast drying. Can be applied more easily and economically than conventional adhesives with watery, thin-flowing consistency.

- Very high adhesion.
- Color Orange: Transparent.

Ref No : 9001 -15 ml Bottle

TRAY

ADHESIVE



A paint on tray adhesive for use with all alginates; it helps to keep the alginate firm in the tray during removal from the mouth. Comes with a convenient brush inside the cap.

Ref No : 5006 -30ml bottle

reviver

SILICON DENTURE RELINING MATERIAL



Long-term soft relining material for chair-side use, A-silicone based, cold curing, biocompatible. The material is directly applied from the cartridge onto the denture base. Strong bonding to the denture acrylic. Also suitable for application with the indirect technique. Primer with colour indicator for better application check. Colour: Gingiva

Ref No : 9006 - 50 ml Cartridge, 7 Mixing Tips, 2x7ml Polish,
5ml Primer

**RELINING
MATERIAL**



www.medicptdental.co.uk



**BITE
REGISTRATION**



**TISSUE
MANAGEMENT**



**TEMPORARY
CEMENTS**

**TRAY
ADHESIVE**



Medicept
DENTAL

MEDICEPT UK LTD.

2nd Floor, Hygeia House, 66 College Road,
Harrow, Middlesex, HA1 1BE, UK
Tel: +44 (0) 1923 251327
Email: info@mediceptdental.co.uk
Website: www.medicptdental.co.uk